


Cyclosporine Level Trough

Order Name: **CYCLO TROU**
Test Number: 2015525
Revision Date: 06/15/2018

TEST NAME	METHODOLOGY	LOINC CODE
Cyclosporine Level Trough	<u>Chemiluminescent Micropartical Immunoassay</u>	53828-0

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	2 mL (1.0)	Whole Blood	EDTA (Lavender Top)	Refrigerated
Instructions	Trough: draw before next dose. Stability: Ambient 24 hours. Refrigerated 7 days. Frozen > 7 days.			

GENERAL INFORMATION	
Testing Schedule	Twice Daily
Expected TAT	1-2 days
Clinical Use	This test is used to monitor blood concentration of Cyclosporine A.
CPT Code(s)	80158
Service Provided By	 labcorp Oklahoma, Inc.