


Creatinine, Serous Fluid

Order Name: **SRS CREAT**
Test Number: 2015850
Revision Date: 06/11/2003

TEST NAME	METHODOLOGY	LOINC CODE
Creatinine, Serous Fluid	<u>Kinetic Alkaline Picrate (Jaffe)</u>	12190-5

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	Serous Fluid	Sterile Screwtop Container	Refrigerated
Instructions	Venous blood is often drawn simultaneously. Note fluid type on requisition and container. Specimen stability: Ambient 8 hours. Refrigerated 7 days.			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	1-2 days
Clinical Use	Useful for sensitive and quantitative measurement of renal function.
CPT Code(s)	82570
Service Provided By	 labcorp Oklahoma, Inc.