


Sodium Urine Random

Order Name: **SODIUM R U**
Test Number: 3001550
Revision Date: 06/04/2021

TEST NAME	METHODOLOGY	LOINC CODE
Sodium Urine Random	<u>Ion-Selective Electrode</u>	2955-3

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	Urine, Random	Sterile Urine container	Refrigerated
Instructions	Random urine collection. No preservative. Keep refrigerated. Specimen stability: Ambient 1 day, Refrigerated 45 days, Frozen 1 year			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	1-2 days
Clinical Use	Useful in monitoring metabolic processes, pituitary function, adrenal function, hydration, proper body pH and regulation of appropriate heart and muscle functions.
CPT Code(s)	84300
Service Provided By	 labcorp Oklahoma, Inc.