


**Uric Acid Urine Random**

Order Name: **URIC R U**  
Test Number: 3001770  
Revision Date: 06/04/2021

TEST NAME	METHODOLOGY	LOINC CODE
Uric Acid Urine Random	Uricase .....	3086-6

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	Urine, Random	Sterile Urine container	Refrigerated
Instructions	Patient should be on normal diet, no alcohol consumption during collection. Random urine collection. <b>Keep refrigerated.</b> Specimen stability: Ambient 24 hours(4 days if pH adjusted >8.0 dropwise with NaOH). Refrigerated 7 days.			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	1-2 days
Clinical Use	Useful in diagnosing and monitoring of therapy in gout.
CPT Code(s)	84560
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.