


Protein Urine Random

Order Name: **PROT R U**
Test Number: 3001950
Revision Date: 05/28/2021

TEST NAME	METHODOLOGY	LOINC CODE
Protein Urine Random	Pyrogallol Red	2888-6

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	Urine, Random	Sterile Urine container	Refrigerated
Instructions	Random urine collection. No preservative. Keep refrigerated. Specimen stability: Ambient 24 hours. Refrigerated 7 days. Frozen 1 month.			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	1-2 days
Clinical Use	Useful for detecting proteinuria and evaluation of renal function.
CPT Code(s)	84156
Service Provided By	 labcorp Oklahoma, Inc.