


Testosterone Total, Adult Male

Order Name: **TESTOS TO**
Test Number: 3602650
Revision Date: 06/15/2018

TEST NAME	METHODOLOGY	LOINC CODE
Testosterone Total, Adult Male	<u>Chemiluminescence Assay</u>	2986-8

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	Serum	Clot Activator SST	Refrigerated
Alternate 1	1 mL (0.5)	Serum	Clot Activator (Red Top, No-Gel)	Refrigerated
Instructions	Not Recommended for Females and Children 13 years and under Preferred to separate serum from cells ASAP or within 2 hours of collection. Stability: Ambient 8 hrs. Refrigerated 7 days. Frozen > 7 days			

GENERAL INFORMATION	
Testing Schedule	Mon-Fri
Expected TAT	1-2 days
Clinical Use	Useful for evaluation of males with erectile dysfunction, gynecomastia, osteoporosis, infertility, delayed or precocious puberty and monitoring replacement therapy. Evaluation of women with hirsutism, virilization and oligomenorrhea and infants with ambiguous genitalia and/or virilizing syndromes.
Notes	For Females and Children 13 years and under use the following test code: Testosterone Total by LCMS [3602651]
CPT Code(s)	84403
Service Provided By	 labcorp Oklahoma, Inc.