

Dehydroepiandrosterone Sulfate (DHEA S)

Order Name: DHEA S
Test Number: 2022725
Revision Date: 12/16/2021

TEST NAME			METHODOLOGY	LOINC CODE	
Dehydroepiandrosterone Sulfate (DHEA S)			Chemiluminescence Assay	2191-5	
SPECIMEN REQUIREMENTS					
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment	
Preferred	1 mL (0.5)	Serum	Clot Activator SST	Refrigerated	
Alternate 1	1 mL (0.5)	Plasma	EDTA (Lavender Top)	Refrigerated	
Alternate 2	1 mL (0.5)	Plasma	Sodium Heparin (Green Top, No-Gel)	Refrigerated	
Instructions	Allow specimen to clot fully, ce	Allow specimen to clot fully, centrifuge, remove serum or plasma from cells. Specimen stability: Refrigerated 8 days, Frozen 2 months.			

GENERAL INFORMATION	
Testing Schedule	Mon-Fri
Expected TAT	1-3 days
Clinical Use	Useful for the diagnosis of congenital adrenal hyperplasia and adrenal carcinoma, determining the cause of hirsutism, virilization and polycystic ovary disease.
CPT Code(s)	82627
Service Provided By	labcorp Oklahoma, Inc.