

Tobramycin


Order Name: **TOBRAMYCIN**

Test Number: 4005550

Revision Date: 06/17/2003

TEST NAME	METHODOLOGY	LOINC CODE
Tobramycin	Enzyme Immunoassay	35670-9

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	Plasma	Lithium Heparin PST (Light Green Top)	Refrigerated
Alternate 1	1 mL (0.5)	Serum	Clot Activator SST	Refrigerated
Instructions	Stability: Ambient 8 hours. Refrigerated 7 days.			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	1-2 days
Clinical Use	Useful for optimizing drug dosage and assessing toxicity.
CPT Code(s)	80200
Service Provided By	<div>labcorp Oklahoma, Inc.</div>