

Drug Screen, Blood

Order Name: **DRUG B SCR**
Test Number: 4300050
Revision Date: 05/03/2021

TEST NAME	METHODOLOGY	LOINC CODE
Acetaminophen Screen Serum	CEDIA	3297-9
Barbiturates Screen	CEDIA	20421-4
Benzodiazepines Screen	CEDIA	46976-7
Ethanol Screen	Enzymatic	20470-1
Salicylate Screen Blood	Colorimetric	4023-8
Tricyclics, Screen		4073-3

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	2 mL (1.0)	Plasma	Lithium Heparin PST (Light Green Top)	Refrigerated
Alternate 1	2 mL (1.0)	Serum	Clot Activator SST	Refrigerated
Instructions	Only original un-opened tube accepted. Specimen stability: Ambient 24 hours, Refrigerated 5 days, Frozen 1 month			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	1-2 days
Clinical Use	Used to evaluate drug abuse, toxicity and therapeutic levels.
CPT Code(s)	80307
Service Provided By	 labcorp Oklahoma, Inc.