

Opiates, Screen Urine

Order Name: **OPIATES U**
Test Number: 4301300
Revision Date: 05/28/2021

TEST NAME	METHODOLOGY	LOINC CODE
Opiates, Screen Urine	CEDIA	3879-4

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	30 mL (10 mL)	Urine, Random	Sterile Urine container	Refrigerated
Instructions	Keep refrigerated. Specimen stability: Refrigerated 5 days, Frozen >5 days			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	1-2 days
CPT Code(s)	80307
Service Provided By	 labcorp Oklahoma, Inc.