


**Tricyclics Screen**

Order Name: **TRICYCL SC**  
Test Number: 4301775  
Revision Date: 01/27/2017

TEST NAME	METHODOLOGY	LOINC CODE
Tricyclics Screen	Enzyme Immunoassay	4073-3

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5)	Serum	Clot Activator (Red Top, No-Gel)	Refrigerated
Instructions	Specimen stability: Ambient 8 hours. Refrigerated 7 days.			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	1-2 days
Clinical Use	Useful for detection of drugs frequently found in drug overdose or used with suicidal intent.
CPT Code(s)	80307
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.