


**Apolipoprotein B**

Order Name: **APO B**  
Test Number: 2015050  
Revision Date: 04/26/2021

TEST NAME	METHODOLOGY	LOINC CODE
<b>Apolipoprotein B</b>	<u>Turbidometric</u> .....	1871-3

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
<b>Preferred</b>	<b>1 mL (0.5)</b>	<b>Serum</b>	<b>Clot Activator SST</b>	<b>Refrigerated</b>
<b>Alternate 1</b>	<b>1 mL (0.5)</b>	<b>Plasma</b>	<b>Lithium Heparin PST (Light Green Top)</b>	<b>Refrigerated</b>
<b>Instructions</b>	<b>Patient Preparation: Freshly drawn fasting specimen.</b> Separate serum/plasma from cells ASAP or within 2 hours of collection Unacceptable Conditions: Hemolyzed specimens. Stability After separation from cells: Ambient 8 hours, Refrigerated 7 days, Frozen 2 months			

GENERAL INFORMATION	
<b>Testing Schedule</b>	Mon - Fri
<b>Expected TAT</b>	2-3 Days
<b>CPT Code(s)</b>	82172
<b>Service Provided By</b>	 <b>labcorp</b> Oklahoma, Inc.