

**Sirolimus (Rapamycin)**

Order Name: **SIROLIMUS**  
Test Number: 3658510  
Revision Date: 06/04/2021

TEST NAME	METHODOLOGY	LOINC CODE
Sirolimus (Rapamycin)	<u>Chemiluminescence Assay</u>	29247-4

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	3 mL (2)	Whole Blood	EDTA (Lavender Top)	Refrigerated
Instructions	Ship and store refrigerated. Shipping at ambient temperature (< 25 C) is acceptable. Optimal time to draw specimen: 1/2 to 1 hour before next oral dose at steady state. Stability: Ambient n/a, Refrigerated 7 days. Frozen > 7 days			

GENERAL INFORMATION	
Testing Schedule	Mon - Fri
Expected TAT	1-3 Days
Notes	*{ Note: 2006 CPT Updated.}
CPT Code(s)	80195
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.