


Apolipoprotein A1

Order Name: **APO A 1**
Test Number: 2015000
Revision Date: 04/26/2021

| TEST NAME | METHODOLOGY | LOINC CODE |
|--------------------------|-------------------------------|------------|
| Apolipoprotein A1 | <u>Turbidometric</u> | 1869-7 |

| SPECIMEN REQUIREMENTS | | | | |
|-----------------------|--|---------------|--|-----------------------|
| Specimen | Specimen Volume (mL) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 1 mL (0.5) | Serum | Clot Activator SST | Refrigerated |
| Alternate 1 | 1 mL (0.5) | Plasma | Lithium Heparin PST (Light Green Top) | Refrigerated |
| Instructions | Patient Preparation: Freshly drawn fasting specimen. Separate serum/plasma from cells ASAP or within 2 hours of collection. Unacceptable Conditions: Hemolyzed specimens. Stability: Ambient 8 hours, Refrigerated 7 days, Frozen 2 months | | | |

| GENERAL INFORMATION | |
|----------------------------|--|
| Testing Schedule | Mon - Fri |
| Expected TAT | 2-3 Days |
| CPT Code(s) | 82172 |
| Service Provided By |  labcorp Oklahoma, Inc. |