


**Helicobacter pylori Antigen Detection Stool**

Order Name: **HPYLORI AG**  
Test Number: 3502325  
Revision Date: 09/21/2022

TEST NAME	METHODOLOGY	LOINC CODE
Helicobacter pylori Antigen Detection Stool	Immunoassay (IA)	80373-4

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	5g (1)	Stool, Random	Stool specimen container	Refrigerated
Instructions	<p><b>Patient Preparation:</b> To avoid false negative results, no antibiotics, proton pump inhibitors, or a bismuth preparation should be administered to patient for 14 days prior to specimen collection.</p> <p>Transfer 5g(1g) stool to an unpreserved stool transport vial and Freeze.</p> <p><b>Unacceptable Conditions:</b> Gastric specimens, swabs or tissue. Specimens in media or preservatives.</p> <p><b>Stability:</b> Ambient: 2 hours; Refrigerated: 72 hours; Frozen: 1 month.</p>			

GENERAL INFORMATION	
Testing Schedule	Sun-Sat
Expected TAT	2 Days
CPT Code(s)	87338
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.