

Cryofibrinogen

Order Name: **CRYOFIBRIN**
Test Number: 5221675
Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
Cryofibrinogen	Precipitation	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	3mL (1mL)	Plasma	Sodium Citrate 3.2% (Blue Top)	See Instructions
Instructions	<p>Specimen Type: Blue-top (sodium citrate) tube, lavender-top (EDTA) tube, or gray-top (sodium fluoride) tube</p> <p>Specimen Storage: Maintain specimen at room temperature.</p> <p>Specimen Collection: Specimen must be drawn in a prewarmed tube and kept at 37(degrees)C. Separate plasma from cells immediately and transfer the plasma into a clean transport tube. LABEL TRANSPORT TUBE "PREWARMED."</p> <p>Specimen Stability: Ambient: Not Available, Refrigerated : Not Available, Frozen: Not Available</p>			

GENERAL INFORMATION	
Expected TAT	3 - 5 days later.
Notes	Labcorp Test Code: 080317
CPT Code(s)	82585
Service Provided By	 labcorp Oklahoma, Inc.