


**Ovarian Antibody Screen with reflex to Titer**

Order Name: **OVARIAN AB**  
Test Number: 5005750  
Revision Date: 06/12/2009

TEST NAME	METHODOLOGY	LOINC CODE
Ovarian Antibody Screen with reflex to Titer	Immunofluorescence Assay (IFA)	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.3)	Serum	Clot Activator SST	Refrigerated
Instructions	<p><b>Notes:</b> 0.3 mL (Note: This volume Does not allow for repeat testing.)</p> <p><b>Specimen Type:</b> Serum from a Gel Tube OR Serum from Red Top Tube</p> <p><b>Specimen Storage:</b> Store frozen or refrigerated.</p> <p><b>Specimen Collection:</b> Separate serum from cells within one hour. Transfer to a plastic transport tube before freezing.</p> <p><b>Special Instructions:</b> Testing is performed at Esoterix Inc. ESECF#55</p> <p><b>Specimen Stability:</b> Ambient: 7 days, Refrigerated : 14 days, Frozen: 60 days</p>			

GENERAL INFORMATION	
Expected TAT	7-14 days
Clinical Use	Ovarian Antibody is found in patients with premature ovarian failure, Addison's disease, and polyendocrinopathy syndromes.
CPT Code(s)	86256
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.