


Group A Streptococcus Culture

Order Name: **C STREP A**
Test Number: 6000250
Revision Date: 06/21/2018

TEST NAME	METHODOLOGY	LOINC CODE
Group A Streptococcus Culture	Culture	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	See Instructions	Swab	Copan eSwab - White (Regular size)	Room Temperature
Alternate 1	See Instructions	Swab	Copan eSwab - Green (Mini-tip)	Room Temperature
Alternate 2	See Instructions	Swab	Double Tipped Aerobic Swab (Red)	Room Temperature
Instructions	For throat collections, swab tonsils or affected area thoroughly and avoid the tongue. Rectal and/or perianal swabs may also be submitted for detecting perianal streptococcal disease. Specimen should be setup within 2 days of collection.			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	2 Days
Clinical Use	Detects Group A Streptococcus by culture.
CPT Code(s)	87081
Service Provided By	 labcorp Oklahoma, Inc.