


**Gonococcus Screen**

Order Name: **C GC**  
Test Number: 6000350  
Revision Date: 09/21/2017

TEST NAME	METHODOLOGY	LOINC CODE
Gonococcus Screen	Culture	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	See Instructions	Swab	Copan eSwab - White (Regular size)	Room Temperature
Alternate 1	See Instructions	Swab	Copan eSwab - Green (Mini-tip)	Room Temperature
Instructions	Obtain culture from mucosal surface with aerobic white swab or green cap minitip swab. Transport at Room Temperature. Culture should be setup within 2 days of collection. Specimens must remain at room temperature, this organism is sensitive to cold temperatures. Refrigerated or frozen specimens are unacceptable.			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	3 Days
Clinical Use	Determines Neisseria gonorrhoeae infections
CPT Code(s)	87081
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.