

Pinworm Exam

Order Name: **C PINWORM**Test Number: 6000600
Revision Date: 12/19/2017

TEST NAME			METHODOLOGY	LOINC CODE
Pinworm Exam			Microscopy	
SPECIMEN REQUIRE	MENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	see instructions	Swube paddle	Swube paddle	Room Temperature
Instructions	Taken with SWUBE (Sticky) Pado	lle. Best taken in middle of t	he night when eggs are laid on the rectal area of p	patient.

GENERAL INFORMATION		
Testing Schedule	Mon - Fri	
Expected TAT	2 Days	
Clinical Use	Detects presence of Enterobius vermicularis (pinworm) larvae	
CPT Code(s)	87172	
Service Provided By	Oklahoma, Inc.	