

Urogenital Culture

Order Name: C UROG RTS
Test Number: 6002005
Revision Date: 09/27/2017

TEST NAME			METHODOLOGY	LOINC CODE
Urogenital Culture			Culture	10352-3
SPECIMEN REQU	JIREMENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred		Swab	Copan eSwab - White (Regular size)	Room Temperature
Alternate 1		Swab	Copan eSwab - Green (Mini-tip)	Room Temperature
Instructions	Specimen should be submitted for culture without delay and transported at room temperature. Unacceptable specimens include those arriving >2 days after collection, refrigerated specimens, specimens in viral transport media, and frozen specimens.			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	4 Days
Clinical Use	Identifies urogenital bacterial pathogens
CPT Code(s)	87070
Service Provided By	labcorp Oklahoma, Inc.