


**KOH Prep for Fungus**

Order Name: **C KOH PR**  
Test Number: 6300100  
Revision Date: 02/21/2020

TEST NAME	METHODOLOGY	LOINC CODE
KOH Prep for Fungus	<a href="#">Microscopy</a>	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	See Instructions	See Instructions	See Instructions	Room Temperature
Instructions	Submit tissue, skin scrapings, or nail pieces in a sterile screw-top cup. For vaginal microscopic exam, send vaginal discharge in a saline wet prep tube.			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	1 Day
Clinical Use	Detection of fungal elements in tissue samples; detection of yeast/fungal elements in vaginal samples
CPT Code(s)	87220
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.