



## Alcohol Plasma/Serum

Order Name: **ALCOHOL**

Test Number: 4000800

Revision Date: 04/26/2021

TEST NAME	METHODOLOGY	LOINC CODE
Alcohol Plasma/Serum	Enzymatic	5643-2

### SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	<b>1 mL (0.1 mL)</b>	<b>Plasma</b>	<b>Lithium Heparin PST (Light Green Top)</b>	<b>Refrigerated</b>
Alternate 1	<b>1 mL (0.1 mL)</b>	<b>Serum</b>	<b>Clot Activator SST</b>	<b>Refrigerated</b>
Instructions	<p>Use betadine or phisoxyhex prep sponge. Do not use alcohol. MUST be submitted in unopened original tube as exposure to air will affect result. Stability: Ambient 2 weeks, Refrigerated 6 months, Frozen 6 months</p>			

### GENERAL INFORMATION

Testing Schedule	Daily
Expected TAT	24hrs
CPT Code(s)	G0480
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.