

Aureobasidium pullulans Allergen

Order Name: **AUREO PULL**
Test Number: **5616780**
Revision Date: 02/11/2013

TEST NAME	METHODOLOGY	LOINC CODE
Aureobasidium pullulans Allergen	ImmunoCAP	15964-0

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	0.4 mL (0.2 mL)	Serum	Clot Activator SST	Room Temperature
Instructions	Specimen Type: Red-top tube or gel-barrier tube, Separate Serum from Cells into a screwtop transport container. Stability Requirements: Room temperature 14 days, Refrigerated 14 days, Frozen 3 months. (Freeze/thaw cycles Stable x3)			

GENERAL INFORMATION	
Expected TAT	3-5 days
CPT Code(s)	86003
Service Provided By	 labcorp Oklahoma, Inc.