

Order Name: GARLIC Test Number: 5581375 Revision Date: 10/23/2017

TEST NAME			METHODOLOGY	LOINC CODE	
Garlic Allergen			ImmunoCAP	15724-8	
SPECIMEN REQUIREMENTS					
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment	
Preferred	0.4 mL (0.2 mL)	Serum	Clot Activator SST	Room Temperature	
Instructions	Specimen Type : Red-top tube or gel-barrier tube, Separate Serum from Cells into a screwtop transport container. Stability Requirements : Room temperature 14 days, Refrigerated 14 days, Frozen 3 months. (Freeze/thaw cycles Stable x3)				

GENERAL INFORMATION	
Expected TAT	3-5 days
CPT Code(s)	86003
Service Provided By	Oklahoma, Inc.