

Penicillin V Allergen

Order Name: PEN V
Test Number: 5510575
Revision Date: 10/01/2023

TEST NAME			METHODOLOGY	LOINC CODE
Penicillin V Allergen			ImmunoCAP	6210-9
SPECIMEN REQUIRE	MENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	0.4 mL (0.2 mL)	Serum	Clot Activator SST	Room Temperature
Instructions	Specimen Type: Red-top tube or gel-barrier tube, Separate Serum from Cells into a screwtop transport container. Stability Requirements: Room temperature 14 days, Refrigerated 14 days, Frozen 3 months. (Freeze/thaw cycles Stable x3)			

GENERAL INFORMATION	
Expected TAT	3-5 days
CPT Code(s)	86003
Service Provided By	Oklahoma, Inc.