

Penicillin Chrysogenium Allergen

Order Name: **PEN CHRYSO**Test Number: 5622400
Revision Date: 12/07/2016

TEST NAME			METHODOLOGY	LOINC CODE
Penicillin Chrysogenium Allergen			ImmunoCAP	15924-4
SPECIMEN REQU	IREMENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	0.4 mL (0.2 mL)	Serum	Clot Activator SST	Room Temperature
Instructions	Specimen Type: Red-top tube or gel-barrier tube,			
Separate Serum from Cells into a screwtop transport container.  Stability Requirements: Room temperature 14 days, Refrigerated 14 days, Frozen 3 months. (Freeze/thaw cycles Stable x				/(h  OtableO)
	Stability Requirements: Roo	gerated 14 days, Frozen 3 months. (Freeze/	rtnaw cycles Stable x3)	

GENERAL INFORMATION	
Expected TAT	3-5 days
CPT Code(s)	86003
Service Provided By	labcorp Oklahoma, Inc.