

St. John Health System Lab Catalog

**Animal and Feather Allergen Panel** 

Order Name: A ANIML PN Test Number: 5616750 Revision Date: 08/01/2023

TEST NAME	METHODOLOGY	LOINC CODE
Cat Dander Allergen	ImmunoCAP	15609-1
Dog Dander Allergen	ImmunoCAP	15685-1
Cow Dander Allergen	ImmunoCAP	15661-2
Horse Dander Allergen	ImmunoCAP	15784-2
Feathers Allergen Mix Panel	ImmunoCAP	31161-3

SPECIMEN REQUIREMENTS					
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment	
Preferred	2 mL (1 mL)	Serum	Clot Activator SST	Room Temperature	
Instructions	<b>Specimen Type</b> : Red-top tube or gel-barrier tube, Separate Serum from Cells into a screwtop transport container. <b>Stability Requirements</b> : Room temperature 14 days, Refrigerated 14 days, Frozen 3 months. (Freeze/thaw cycles Stable x3)				

GENERAL INFORMATION		
Expected TAT	3-5 Days	
Notes	Performed at Labcorp Phoenix as individual allergens	
CPT Code(s)	86003x5	
Service Provided By	O labcorp	