


Animal and Feather Allergen Panel

Order Name: **A ANIML PN**
Test Number: **5616750**
Revision Date: **08/01/2023**

| TEST NAME | METHODOLOGY | LOINC CODE |
|-----------------------------|-------------|------------|
| Cat Dander Allergen | ImmunoCAP | 15609-1 |
| Dog Dander Allergen | ImmunoCAP | 15685-1 |
| Cow Dander Allergen | ImmunoCAP | 15661-2 |
| Horse Dander Allergen | ImmunoCAP | 15784-2 |
| Feathers Allergen Mix Panel | ImmunoCAP | 31161-3 |

| SPECIMEN REQUIREMENTS | | | | |
|-----------------------|---|---------------|--------------------|-----------------------|
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 2 mL (1 mL) | Serum | Clot Activator SST | Room Temperature |
| Instructions | Specimen Type: Red-top tube or gel-barrier tube, Separate Serum from Cells into a screwtop transport container. Stability Requirements: Room temperature 14 days, Refrigerated 14 days, Frozen 3 months. (Freeze/thaw cycles Stable x3) | | | |

| GENERAL INFORMATION | |
|---------------------|--|
| Expected TAT | 3-5 Days |
| Notes | Performed at Labcorp Phoenix as individual allergens |
| CPT Code(s) | 86003x5 |
| Service Provided By |  labcorp Oklahoma, Inc. |