

Burkitt's Lymphoma/NHL/ALL, IGH/MYC, t(8;14) by FISH

Order Name: **BURKI FSH**Test Number: 9117900
Revision Date: 10/01/2022

TEST NAME			METHODOLOGY	LOINC CODE
Burkitt's Lymphoma/NHL/ALL, IGH/MYC, t(8;14) by FISH			Fluorescence in Situ Hybridization	
SPECIMEN REQU	JIREMENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	5 mL (1 mL)	Bone Marrow	Sodium Heparin (Green Top, No-Gel)	Room Temperature
Alternate 1	5 mL (3 mL)	Whole Blood	Sodium Heparin (Green Top, No-Gel)	Room Temperature
Instructions	Send specimen ASAP, Keep a	Send specimen ASAP, Keep at room temperature! (DO NOT FREEZE). Frozen samples will be rejected.		

GENERAL INFORMATION	
Testing Schedule	Mon-Fri
Expected TAT	3-5 Days
Clinical Use	Useful to diagnose Burkitt's type lymphoma and some cases of acute lymphoblastic leukemia.
CPT Code(s)	88377
Service Provided By	labcorp Oklahoma, Inc.