

limit of quantitation for the assay.

Oklahoma, Inc.

Labcorp Test Code: 605059

86001

Notes

CPT Code(s)

Service Provided By

Order Name: COCONUT G Test Number: 5500471 Revision Date: 12/14/2022

TEST NAME			METHODOLOGY	LOINC CODE
Coconut IgG	G Enzyme immunoassay (EIA)			
SPECIMEN REQUI	REMENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	0.3 mL (0.1 mL)	Serum	Clot Activator SST	Room Temperature
Instructions	Specimen Type: Red-top tube or gel-barrier tube			
	Specimen Storage: Room temperature			
	Specimen Collection: Not Available			
	Special Instructions: Values obtained with different assay methods should not be used interchangeably in serial testing. It is recommended that only			
	one assay method be used consistently to monitor each patient's course of therapy.			
	Specimen Stability: Ambient: Not Available, Refrigerated : Not Available, Frozen: Not Available			
GENERAL INFORM	IATION			
Expected TAT	2 - 4 days			
Clinical Use	The clinical utility of food-specific IgG tests has not been established. These tests can be used in special clinical situations to select foods for			

evaluation by diet elimination and challenge in patients who have food-related complaints. It should be recognized that the presence of foodspecific IgG alone cannot be taken as evidence of food allergy and only indicates immunologic sensitization by the food allergen in question. This test should only be ordered by physicians who recognize the limitations of the test. The reference range listed on the report is the lower