Stemphylium herbarum IgG

## St. John Health System Lab Catalog

Order Name: STEMPHYL G Test Number: 5500449

Revision Date: 12/14/2022

TEST NAME		М	ETHODOLOGY	LOINC CODE	
Stemphylium herbarum IgG		Ë	Enzyme immunoassay (EIA)	35553-7	
SPECIMEN REQUIREMENTS					
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment	
Preferred	0.3 mL (0.1 mL)	Serum	Clot Activator (SST or Red No-Gel)	Room Temperature	
Instructions	Specimen Type: Red-top tube or gel-barrier tube				
	Specimen Storage: Room temperature  Specimen Collection: Not Available  Special Instructions: Values obtained with different assay methods should not be used interchangeably in serial testing. It is recommended that only				
	one assay method be used consistently to monitor each patient's course of therapy.				
	Specimen Stability: Ambient: Not Available, Refrigerated : Not Available, Frozen: Not Available				

GENERAL INFORMATION		
Expected TAT	2 - 4 days	
Clinical Use	M010-lgG Stemphylium herbarum  Although there have been many publications concerning the measurement of allergen-specific lgG, the clinical utility of such tests has not been established except in special situations. Thus, the quantitative lgG test should only be ordered by specialists who recognize the limitations of the test. The normal reference ranges reported represent the expected results for individuals who have no unusual exposure and have not been immunized with the indicated allergen. The ranges reported have no disease-associated significance.	
CPT Code(s)	86001	
Service Provided By	labcorp Oklahoma, Inc.	