St. John Health System Lab Catalog

Penicillium Chrysogenym/notatum IgG

Order Name: PEN CHRY G
Test Number: 5500433
Revision Date: 12/14/2022

TEST NAME		METHO	DDOLOGY	LOINC CODE	
Penicillium Chrysogenym/notatum IgG		Enzym	ne immunoassay (EIA)	26957-	
SPECIMEN REQUIREMENTS					
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment	
Preferred	0.3 mL (0.1 mL)	Serum	Clot Activator SST	Room Temperature	
Instructions	Specimen Type: Red-top tube or gel-barrier tube Specimen Storage: Room temperature Specimen Collection: Not Available Special Instructions: Values obtained with different assay methods should not be used interchangeably in serial testing. It is recommended that only one assay method be used consistently to monitor each patient's course of therapy. Specimen Stability: Ambient: Not Available, Refrigerated: Not Available, Frozen: Not Available				

GENERAL INFORMATION	
Expected TAT	2 - 4 days
Notes	Labcorp Test Code: 605213
CPT Code(s)	86001
Service Provided By	labcorp Oklahoma, Inc.