

Acid Phosphate Total

Order Name: ACID PHS T
Test Number: 2000280
Revision Date: 12/12/2022

TEST NAME		MET	THODOLOGY	LOINC CODE	
Acid Phosphate Total		Enz	zymatic	12173-1	
SPECIMEN REQUIREMENTS					
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment	
Preferred	1.5 mL (0.8 mL)	Serum	Clot Activator (Red Top, No-Gel)	Frozen	
Instructions	Allow specimen to clot completely. Transfer 1.5mL of serum into an aliquot tube and freeze. CRITICAL FROZEN. Separate specimens must be submitted when multiple tests are ordered. Unacceptable Conditions: Plasma. Non-frozen specimens. Hemolyzed specimens. Stability After separation from cells: Ambient: Unacceptable; Refrigerated: Unacceptable; Frozen: 1 month				

GENERAL INFORMATION		
Expected TAT	5 - 8 days	
Clinical Use	This test is used in the evaluation of prostate cancer, Padget's disease, hyperthyroidism, and Gaucher's disease.	
Notes	Labcorp Test Code: 829642	
CPT Code(s)	84060	
Service Provided By	labcorp Oklahoma, Inc.	