


**Galectin-3**

Order Name: **GALECTIN 3**  
Test Number: 2015057  
Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
<b>Galectin-3</b>	<u>Quantitative Enzyme Immunoassay</u>	62419-7

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
<b>Preferred</b>	<b>0.8 mL (0.4)</b>	<b>Plasma</b>	<b>EDTA (Lavender Top)</b>	<b>Refrigerated</b>
<b>Alternate 1</b>	<b>0.8 mL (0.4)</b>	<b>Serum</b>	<b>Clot Activator (Red Top, No-Gel)</b>	<b>Refrigerated</b>
<b>Instructions</b>	<b>Notes:</b> 0.4 mL (Note: This volume Does NOT allow for repeat testing.) <b>Specimen Type:</b> Red-top tube, gel-barrier tube OR lavender-top (EDTA) tube <b>Specimen Storage:</b> Room temperature Stability: Room temperature up to 22 days Refrigerated up to 22 days Frozen up to 12 months Stable x9 <b>Specimen Collection:</b> Transfer separated serum or plasma to a plastic transport tube. <b>Specimen Stability:</b> Ambient: Not Available, Refrigerated : Not Available, Frozen: Not Available			

GENERAL INFORMATION	
<b>Expected TAT</b>	4 - 6 days
<b>Notes</b>	Labcorp Test Code: 004110
<b>CPT Code(s)</b>	82777
<b>Service Provided By</b>	 <b>labcorp</b> Oklahoma, Inc.