

Watermelon IgG

Instructions

Order Name: WATERMLN G

Test Number: 5500649
Revision Date: 12/14/2022

| TEST NAME | | N | IETHODOLOGY | LOINC CODE |
|------------------|-----------------------|---------------|--------------------|-----------------------|
| Watermelon IgG | | | | |
| SPECIMEN REQUIRE | MENTS | | | |
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 0.3 mL (0.1 mL) | Serum | Clot Activator SST | Room Temperature |
| | | | | |

Specimen Type: Red-top tube or gel-barrier tube

Specimen Storage: Room temperature **Specimen Collection:** Not Available

Special Instructions: Values obtained with different assay methods should not be used interchangeably in serial testing. It is recommended that only

one assay method be used consistently to monitor each patient's course of therapy.

Specimen Stability: Ambient: Not Available, Refrigerated : Not Available, Frozen: Not Available

| GENERAL INFORMATION | | |
|---------------------|---------------------------|--|
| Expected TAT | 2 - 4 days | |
| Notes | Labcorp Test Code: 605054 | |
| CPT Code(s) | 86001 | |
| Service Provided By | labcorp Oklahoma, Inc. | |