


**Coumadin (Warfarin) Level**

Order Name: **Coumadin**  
Test Number: 3606450  
Revision Date: 08/31/2023

TEST NAME	METHODOLOGY	LOINC CODE
Coumadin (Warfarin) Level	<u>HPLC/Tandem Mass Spectrometry</u>	4098-0

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.3 mL)	Serum	Clot Activator (Red Top, No-Gel)	Room Temperature
Alternate 1	1 mL (0.3 mL)	Plasma	EDTA (Lavender Top)	Room Temperature
Instructions	<p><b>Specimen Type:</b> Red-top tube, lavender-top (EDTA) tube <b>DO NOT USE A GEL-BARRIER TUBE.</b></p> <p>The use of gel-barrier tubes is not recommended due to slow absorption of the drug by the gel.</p> <p><b>Specimen Collection:</b> Promptly centrifuge and separate serum or plasma into a plastic screw capped vial using approved guidelines.</p> <p><b>Specimen Stability:</b> Ambient: 30 days, Refrigerated : 30 days, Frozen: 15 months</p>			

GENERAL INFORMATION	
Expected TAT	3 - 4 days
CPT Code(s)	80375
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.