


## Butalbital Level (Fiorinal)

Order Name: **Butalbital**  
Test Number: 3608250  
Revision Date: 11/13/2024

TEST NAME	METHODOLOGY	LOINC CODE
Butalbital Level (Fiorinal)	<u>Liquid Chromatography/Tandem Mass Spectrometry</u>	6895-7

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	2.0mL (0.5mL)	Serum	Clot Activator (Red Top, No-Gel)	Refrigerated
Alternate 1	2.0mL (0.5mL)	Plasma	Sodium Heparin (Green Top, No-Gel)	Refrigerated
Instructions	<p><b>Specimen:</b> 2mL (0.5mL) Serum or Plasma - Refrigerated</p> <p><b>Container:</b> Red-Top (no additive) Tube or Green-Top (heparin) Tube. <b>Gel- barrier tubes are not recommended.</b></p> <p>Testing referred to MedTox Laboratories Inc. MX#22</p> <p><b>Collection:</b> Transfer separated serum or plasma to a plastic transport tube.</p> <p><b>Storage Instructions:</b> Submission/transport (&lt;3 days) Room temperature. For storage beyond three days, specimen should be Refrigerated or Frozen.</p>			

GENERAL INFORMATION	
Expected TAT	Set up Monday thru Friday on 1st shift; reported 1-2 Days from setup
Clinical Use	Also known as: Axocet, Lotunate, and Sandoptal
CPT Code(s)	80345
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.