

Order Name: GIARDIA AB Test Number: 3805575 Revision Date: 12/12/2022

| TEST NAME | | | METHODOLOGY | LOINC CODE |
|-----------------------|--|---------------|----------------------------------|-----------------------|
| Giardia Antibody IgG | | | Radial Immunodiffusion | |
| SPECIMEN REQUIREMENTS | | | | |
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 1.0mL (0.5mL) | Serum | Clot Activator SST | Frozen |
| Alternate 1 | 1.0mL (0.5mL) | Serum | Clot Activator (Red Top, No-Gel) | Frozen |
| Instructions | Specimen Type: Red-top tube or gel-barrier tube Specimen Storage: FREEZE Stability: Refrigerate up to 7 days Frozen up to 60 days Specimen Collection: Serum should be separated from cells within 2 hours of venipuncture. Send serum in a plastic transport tube. Special Instructions: Testing referred to BioAgilytix Diagnostics Lab E=#- Specimen Stability: Ambient: Not Available, Refrigerated : Not Available, Frozen: Not Available | | | |

| GENERAL INFORMATION | | |
|---------------------|---------------------------|--|
| Expected TAT | 11-14 days | |
| Notes | Labcorp Test Code: 800602 | |
| CPT Code(s) | 86674 | |
| Service Provided By | O labcorp | |