

## Hepatitis B Virus Genotype by Sequencing

Order Name: HBV Genotyp Seq

Test Number: 3488232 Revision Date: 02/13/2017

| TEST NAME                                |                                                                                                                                               |                                          | METHODOLOGY                 | LOINC CODE            |  |
|------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------|-----------------------------|-----------------------|--|
| Hepatitis B Virus Genotype by Sequencing |                                                                                                                                               |                                          | Polymerase Chain Reaction   |                       |  |
| SDECIMEN DECLUDES                        | MENTS                                                                                                                                         |                                          |                             |                       |  |
| SPECIMEN REQUIREM                        | WENTS                                                                                                                                         |                                          |                             |                       |  |
| Specimen                                 | Specimen Volume (min)                                                                                                                         | Specimen Type                            | Specimen Container          | Transport Environment |  |
| Preferred                                | 1.2 mL (600 uL)                                                                                                                               | Serum                                    | Clot Activator SST          | Frozen                |  |
| Instructions                             | Specimen Type: Gel-barrier tube, yellow-top (ACD) tube, lavender-top (EDTA) tube or PPT.                                                      |                                          |                             |                       |  |
|                                          | Specimen Storage: FREEZE                                                                                                                      |                                          |                             |                       |  |
|                                          | Specimen Collection: Centrifuge sample within six hours of collection, transfer plasma/serum to a polypropylene screw-cap tube and FREEZE. To |                                          |                             |                       |  |
|                                          | avoid delays in turnaround time,                                                                                                              | sts on frozen samples, PLEASE SUBMIT SEF | PARATE FROZEN SPECIMENS FOR |                       |  |
|                                          | EACH TEST REQUESTED.                                                                                                                          |                                          |                             |                       |  |

| GENERAL INFORMATION |                           |  |  |
|---------------------|---------------------------|--|--|
| Expected TAT        | 7-10 days                 |  |  |
| Notes               | Labcorp Test Code: 551710 |  |  |
| CPT Code(s)         | 87912                     |  |  |
| Service Provided By | labcorp Oklahoma, Inc.    |  |  |