


## EBV (Epstein Barr Virus) Panel

Order Name: **EBV Pnl**  
Test Number: 5580201  
Revision Date: 04/11/2023

TEST NAME	METHODOLOGY	LOINC CODE
EBV Early Antigen IgG	<a href="#">Chemiluminescence Assay</a>	24007-7
EBV Early Antigen IgG Interpretation		40752-8
EBV VCA IgG	<a href="#">Chemiluminescence Assay</a>	5157-3
EBV VCA IgG Interpretation		24114-1
EBV VCA IgM	<a href="#">Chemiluminescence Assay</a>	5159-9
EBV VCA IgM Interpretation		24115-8
EBV Nuclear IgG	<a href="#">Chemiluminescence Assay</a>	30083-0
EBV Nuclear IgG Interpretation		5156-5

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	3 mL (0.75mL)	Serum	Clot Activator SST	Refrigerated or Frozen
Instructions	Specimen Stability: Room Temperature: Unacceptable, Refrigerated: 7 days, Frozen 30 Days			

GENERAL INFORMATION	
Testing Schedule	Mon-Sat
Expected TAT	1-2 Days
Notes	EB Virus IgM IFA confirmation testing will no longer be performed due to supply chain availability. The results for this test are not, by themselves, diagnostic and should be considered in association with clinical presentation and clinical history.
CPT Code(s)	86663, 86664, 86665x2
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.