

**Alpha-Galactosidase, Serum (Fabry Disease)**

Order Name: **Alpha-Galactos**  
Test Number: 6904743  
Revision Date: 12/10/2022

TEST NAME	METHODOLOGY	LOINC CODE
Alpha-Galactosidase, Serum (Fabry Disease)	Quantitative Fluorometry	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	2 mL (0.2 mL)	Serum	Clot Activator SST	Frozen
Instructions	<p><b>Specimen Type:</b> Red-top tube or gel-barrier tube <b>Specimen Storage:</b> FREEZE <b>Specimen Collection:</b> Serum must be separated from cells within 45 minutes of venipuncture. Send serum in a plastic transfer tube. <b>Special Instructions:</b> Testing referred to Mayo Medical Laboratory ML#AGAS Was 8784 <b>Specimen Stability:</b> Ambient: Not Available, Refrigerated : 24 hours, Frozen: 14 days</p>			

GENERAL INFORMATION	
Expected TAT	7-14 days
Notes	Labcorp Test Code: 080291
CPT Code(s)	82657
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.