


Deoxypyridinoline Crosslinks (Pyrilinks-D)

Order Name: **Pyrilinks-D**
Test Number: 6904991
Revision Date: 01/28/2025

| TEST NAME | METHODOLOGY | LOINC CODE |
|--|---------------------------------------|------------|
| Deoxypyridinoline Crosslinks (Pyrilinks-D) | <u>Immunochemiluminometric (ICMA)</u> | |

| SPECIMEN REQUIREMENTS | | | | |
|-----------------------|---|-------------------|-------------------------|-----------------------|
| Specimen | Specimen Volume (mL) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 10 mL (1 mL) | Urine, First Void | Sterile Urine container | Refrigerated |
| Instructions | <p>Specimen Type: Plastic urine container Specimen Storage: Refrigerated Specimen Collection: Collection before 10 AM is recommended (eg, first morning void). Do NOT add preservative. Special Instructions: The account must submit the patient's Social Security number to monitor. Specimen Stability: Ambient: 14 days, Refrigerated : 14 days, Frozen: 14 days</p> | | | |

| GENERAL INFORMATION | |
|---------------------|--|
| Expected TAT | 4-5 days |
| CPT Code(s) | 82570, 82523 |
| Service Provided By |  labcorp Oklahoma, Inc. |