


HLA-B*5701 Typing (Abacavir hypersensitivity)

Order Name: **HLA-B5701**
Test Number: 6905039
Revision Date: 12/09/2022

TEST NAME	METHODOLOGY	LOINC CODE
HLA-B*5701 Typing (Abacavir hypersensitivity)	<u>Polymerase Chain Reaction</u>	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	7 mL (3 mL)	Whole Blood	EDTA (Lavender Top)	Room Temperature
Instructions	<p>Notes: 3 mL whole blood or four (4) buccal swabs (Note: This volume Does NOT allow for repeat testing.)</p> <p>Specimen Type: Lavender-top (EDTA) tube or four (4) buccal swabs in a sealed envelope (buccal swab kit). When submitting buccal swabs please use a buccal swab kit provided by LabCorp. To obtain the buccal kit, please telephone 800-533-1037.</p> <p>Container Detail: WhlBld</p> <p>Specimen Storage: Maintain blood specimen at room temperature. Keep buccal swabs dry and at room temperature.</p> <p>Specimen Collection: Not Available</p> <p>Special Instructions: If you have questions, please telephone 800-533-1037 (HLA customer service) for assistance selecting the proper HLA test for the patient.</p> <p>Specimen Stability: Ambient: Not Available, Refrigerated : Not Available, Frozen: Not Available</p>			

GENERAL INFORMATION	
Expected TAT	3 - 7 days
Notes	Labcorp Test Code: 006926
CPT Code(s)	81381
Service Provided By	 labcorp Oklahoma, Inc.