

Trichomonas Urine - TMA

Order Name: **Trich Urine**
Test Number: **6910123**
Revision Date: **07/11/2017**

TEST NAME	METHODOLOGY	LOINC CODE
Trichomonas Urine	Transcription-Mediated Amplification	70167-2

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	2 mL	Urine, First Void	APTIMA Urine Transport Tube	Refrigerated
Instructions	First catch urine. Patient should not have urinated for at least 1 hr prior to collection. Female patients should not cleanse the labial area prior to collection. Transfer 2 mL urine to APTIMA Urine Specimen Transport Tube within 24hrs of collection. Liquid level must be between fill lines on tube. Stability: Refrigerated = 30 days.			
Please contact your account rep if training and APTIMA Urine Specimen Transport supplies have not been supplied to your practice.				

GENERAL INFORMATION	
Testing Schedule	Mon-Fri
Expected TAT	1-2 days after set-up
CPT Code(s)	87661
Service Provided By	 labcorp Oklahoma, Inc.