

Volatiles, Whole Blood

Order Name: Volatiles WB
Test Number: 4302200
Revision Date: 07/07/2025

TEST NAME	METHODOLOGY	LOINC CODE
Acetone Level	Gas chromatography (GC)	5568-1
Methanol Level	Gas chromatography (GC)	5693-7
Isopropanol Level	Gas chromatography (GC)	5669-7
Ethanol Level	Gas chromatography (GC)	56478-1

SPECIMEN REQUIREMENTS					
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment	
Preferred	7mL (0.5mL)	Whole Blood	Sodium Floride (Gray)	Room Temperature	
Alternate 1	7mL (0.5mL)	Serum	No Additive Clot (Red Top, No-Gel, Plastic)	Room Temperature	
Instructions	Specimen: 7 mL (0.5 mL) Whole blood or serum - Room Temperature Container: Gray-top (sodium fluoride) tube (preferred) or red-top tube. Submit original unopened tube. Stability Requirements: Room temperature 14 days, Refrigerated 14 days, Frozen 14 days (Freeze/thaw cycles Stable x3) Special Instructions: This test number is not intended for workplace testing and does not comply with state regulatory workplace testing programs. Do not prepare venipuncture site with alcohol or remove stopper from tube.				

GENERAL INFORMATION		
Expected TAT	2-3 days	
CPT Code(s)	80320	
Service Provided By	labcorp Oklahoma, Inc.	