

Cedar/Red Tree IgE

Order Name: Red Cedar Tree

Test Number: 5613499
Revision Date: 12/14/2022

| TEST NAME             |  |               | METHODOLOGY        | LOINC CODE            |  |
|-----------------------|--|---------------|--------------------|-----------------------|--|
| Cedar/Red Tree IgE    |  |               | ImmunoCAP          |                       |  |
| SPECIMEN REQUIREMENTS |  |               |                    |                       |  |
| Specimen              | Specimen Volume (min)  | Specimen Type | Specimen Container | Transport Environment |  |
| Preferred             | 0.25 mL (0.1 mL)   | Serum         | Clot Activator SST | Room Temperature      |  |
| Instructions          | Specimen Type: Red-top tube or gel-barrier tube Specimen Storage: Room temperature Specimen Collection: Not Available Specimen Stability: Ambient: 14 days, Refrigerated : 14 days, Frozen: 3 months |               |                    |                       |  |

| GENERAL INFORMATION |                           |  |
|---------------------|---------------------------|--|
| Expected TAT        | 2 - 4 days                |  |
| Notes               | Labcorp Test Code: 605616 |  |
| CPT Code(s)         | 86003                     |  |
| Service Provided By | labcorp Oklahoma, Inc.    |  |