

Biotinidase, Serum (with Paired Normal Control)

Order Name: **Biotinidase**
Test Number: **5613879**
Revision Date: **12/12/2022**

TEST NAME	METHODOLOGY	LOINC CODE
Biotinidase, Serum (with Paired Normal Control)	Spectrophotometry	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (1 mL)	Serum	Clot Activator (Red Top, No-Gel)	Frozen
Instructions	<p>Specimen Type: Red-top tube Specimen Storage: FREEZE Specimen Collection: Separate serum from cells within 30 minutes of venipuncture. Transfer specimen to plastic transport tube before freezing. To avoid delays in turnaround time when requesting multiple tests on frozen samples, please submit separate frozen specimens for each test requested. Special Instructions: ***NOTE - DO NOT USE THIS TEST CODE FOR BIOTIN/VIT B7*** TESTING Do not modify test inquiry names without permission of CMBP.</p>			

GENERAL INFORMATION	
Expected TAT	4-10 days
Notes	Labcorp Test Code: 402362
CPT Code(s)	82261
Service Provided By	 labcorp Oklahoma, Inc.