


**Fentanyl Urine Screen with Reflex to Quant (Send Out)**

Order Name: **Fentanyl Ur Scr**  
Test Number: **5613589**  
Revision Date: 04/15/2021

TEST NAME	METHODOLOGY	LOINC CODE
Fentanyl Urine Screen with Reflex to Quant (Send Out)	Immunoassay (IA)	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	20 mL	Urine, Random	Sterile Urine container	Room Temperature
Instructions	<p><b>Specimen Type:</b> Plastic urine container</p> <p><b>Specimen Storage:</b> Maintain specimen at room temperature. Stability: If arrival at lab will extend beyond 7 days then refrigerate.</p> <p><b>Special Instructions:</b> This assay is intended for pain management. It is not intended for workplace testing and does not comply with state regulatory workplace testing programs.</p>			

GENERAL INFORMATION	
Expected TAT	1-4 days
Clinical Use	<p>If the specimen screens positive, then Confirmation/Quantitation by LC-MS/MS</p> <p><b>Reflex Order code: (5613591) .Fentanyl &amp; Metabolite, Uri, Quant</b></p> <p><b>Result code: (5613591) Norfentanyl Quant</b></p> <p><b>Result code: (5613877) Fentanyl Quant</b></p>
CPT Code(s)	80307
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.