

**Flunitrazepam and Metabolites, Quantitative, Serum**

Order Name: **FLUNTR SCR**  
Test Number: 4013050  
Revision Date: 01/04/2023

TEST NAME	METHODOLOGY	LOINC CODE
Flunitrazepam and Metabolites, Quantitative, Serum	<a href="#">Liquid Chromatography/Tandem Mass Spectrometry</a>	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	5 mL (2 mL)	Serum	Clot Activator (Red Top, No-Gel)	Refrigerated
Instructions	<p><b>Specimen Type:</b> Red-top tube OR green-top (heparin) tube. Gel-barrier tubes are NOT recommended.</p> <p><b>Specimen Storage:</b> Room temperature is adequate for shipping and/or short-term storage (three days); specimens may also be refrigerated or frozen.</p> <p><b>Specimen Collection:</b> Serum or plasma should be separated from cells within two hours of venipuncture. Submit serum or plasma in a plastic transport tube.</p> <p><b>Special Instructions:</b> Testing referred to MedTox Laboratories MX#438</p>			

GENERAL INFORMATION	
Expected TAT	10-14 days
Notes	Labcorp Test Code: 808443
CPT Code(s)	80346 (HCPCS: G0480)
Service Provided By	 <b>labcorp</b> Oklahoma, Inc.